

The **brakes** are located on both sides, above the front wheels. First look carefully which way you have to push the switch to unlock. You cannot push the chair without first unlocking. Please **DO NOT** use YOUR FEET to (un)lock the brakes, you might damage them.

If you need to lift the chair manually, please lift the wheelchair from the **base only**. The chair is very heavy, so lift it by grabbing the inside of the front tires and grabbing the metal bars on the rear tires.

You need 4 people to safely lift the wheelchair, so **one on each tire**. The brakes must be released to be able to push or pull the seat.

Please, consult

(*Phone:*) if there are any problems (un)loading the wheelchair! This chair is of vital importance to us. And can easily be damaged by lifting incorrectly.

Please do not place the seat on a moving loading belt without supervision. Do not push or lift using the handrails or seat section. Separate the seat from all luggage that can fall on it.

Thank you

THIS DOCUMENT IS PROVIDED BY ROLL ON ADVENTURES